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Adaway JE, Keevil BG, Therapeutic drug monitoring and LC-MS/MS, *J. Chromatogr. B Analyt Technol Biomed Life Sci.*883–884:33–49 (2012)

Agrawal YP, Cid M, Westgard S, Parker TS, Jaikaran R, Levine DM. Transplant patient classification and tacrolimus assays: More evidence of the need for assay standardization. *Ther. Drug Monit.* 36:706-709 (2014)

Armbruster DA, Pry T. Limit of blank, limit of detection and limit of quantitation. *Clin Biochem Rev. Suppl* 1:S49-52 (2008)

Aranda JV, Beharry K, Valencia GB, Natarajan G, Davis J. Caffeine impact on neonatal morbidities. *J Matern Fetal Neonatal Med Suppl.* 3: 20-23 (2010)

Annesley TM, McKeown DA, Holt DW, Mussell C, Champamaud E, Harter L, Calton LJ, Maison DS. Standardization of LC-MS for therapeutic drug monitoring of tacrolimus. *Clin. Chem.* 59:1630-1637 (2013)

Apple FS. Clinical and analytical standardization issues confronting cardiac troponin I. *Clin. Chem* 45:18-20 (1999)

Apple FS. Standardization of cardiac markers. *Scand J Clin Lab Invest Suppl.* 240:107-111 (2005)

Apple FS, Jesse RL, Newby LK, Wu AH, Christenson RH, Cannon CP, Francis G, Morrow DA, Ravkilde J, Storrow AB, Tang W; IFCC Committee on Standardization of Markers of Cardiac Damage, Jaffe AS, Mair J, Ordonez-Llanos J, Pagani F, Panteghini M, Tate J; National Academy of Clinical Biochemistry. National Academy of Clinical Biochemistry and IFCC Committee for Standardization of Markers of Cardiac Damage Laboratory Medicine Practice Guidelines: analytical issues for biochemical markers of acute coronary syndromes. *Clin Chem.*53:547-551 (2007)

Aranda JV, Beharry K, Rex RJ, Johannes RJ, Charest-Boule I. Caffeine enzyme immunoassay in neonatal and pediatric drug monitoring. *Ther Drug Monit* 9:97-103 (1987)

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Becker S, Thiery J, Ceglarek U. Evaluation of a novel commercial assay for the determination of cyclosporine A, tacrolimus, sirolimus, and everolimus by Liquid Chromatography-Tandem Mass Spectrometric assay. *Ther. Drug Monit.* 35:129-132 (2013)

Berger P., Laphorn A Standardization of epitopes for human chorionic gonadotropin (hCG) immunoassays. *Curr. Med. Chem.* 23:3481-3494 (2016)

BIO-RAD EQAS PROGRAMME USER GUIDE (<http://www.qcnet.com>)

Budd JR, Weykamp C, Raj R, Mackenzie F, Ceriotti F, Greenberg N, Camara JE, Schimmel H, Vesper HW, Keller T, Delatour V, Panteghini M, Burns C, Miller G for the Working Group on Commutability. IFCC

Working Group Recommendations for assessing commutability Part 3: Using the calibration effectiveness of a reference material. *Clin. Chem.* 64:3 (2018)

Buoro S, Da Rin G, Fanelli A, Lippi G. Harmonization of interpretative comments in laboratory hematology reporting: the recommendations of Working Group on Diagnostic Hematology of the Italian Society of Clinical Chemistry and Clinical Molecular Biology (WGDH-SIBioC). *Clin Chem Lab Med.* 2018 pii: /j/cclm.ahead-of-print/cclm-2017-0972/cclm-2017-0972.xml.

Brereton M, McCafferty R, Marsden K, Kawai Y, Etzell J, Ermens A; international council for standardization in haematology. Recommendation for standardization of haematology reporting units used in the extended blood count. *Int J Lab Hematol.* 38:472-482 (2016)

Cardoso E, Mercier T, Wagner AD, Homicsko K, Michielin O, Ellefsen-Lavoie K, Cagnon L, Diezi M, Buclin T, Widmer N, Csajka C, Decosterd L. Quantification of the next-generation oral anti-tumor drugs dabrafenib, trametinib, vemurafenib, cobimetinib, pazopanib, regorafenib and two metabolites in human plasma by liquid chromatography-tandem mass spectrometry. *Chromatogr B Analyt Technol Biomed Life Sci.* 1083:124-136 (2018)

Carobene A, Guerra E, Ceriotti F. A mechanism-based way to evaluate commutability of control materials for enzymatic measurements. The example of gamma-glutamyltransferase. *Clin Chim Acta.* 424:153-8 (2013)  
Ceriotti F., Hinzmann R, Panteghini M. Reference intervals: the way forward. *Ann. Clin. Biochem.* 46:8-17 (2006)

Cattozzo G, Guerra E, Ceriotti F, Franzini C; Enzyme Working Group of the Italian Society of Clinical Biochemistry and Clinical Molecular Biology (SIBioC). Commutable calibrator with value assigned by the IFCC reference procedure to harmonize serum lactate dehydrogenase activity results measured by 2 different methods. *Clin Chem.* 54:1349-1355 (2008)

Chacko S, Haseeb S, Glover BM, Wallbridge D, Harper A. The role of biomarkers in the diagnosis of acute coronary syndrome. *Future Sci. OA* 4(1), FS0251 (2018)

Chan EK, Damoiseaux J, Carballo OG, Conrad K, de Melo Cruvinel W, Francescantonio PL, Fritzler MJ, Garcia-De La Torre I, Herold M, Mimori T, Satoh M, von Mühlen CA, Andrade LE. Report of the First International Consensus on Standardized Nomenclature of Antinuclear Antibody HEp-2 Cell Patterns 2014-2015. *Front Immunol.* 6:412 (2015)

Christenson RH, Vaidya H, Landt Y, Bauer RS, Green SF, Apple FA, Jacob A, Magnuson GR, Nag S, Wu AH, Azzazy HM. Standardization of creatine kinase –MB (CK-MB) mass assays: the use of recombinant CK-MB as a reference material. *Clin.Chem.* 45: 1414-143 (1999)

Christenson RH, Duh SH, Apple FS, Bodor G, Bunk DM, Dalluge J, Panteghini M, Potter JD, Welch MJ, Wu AH, Kahn SE. Standardization of cardiac troponin I assays: round robin of ten candidate reference materials. *Clin. Chem.* 47: 431-437 (2001)

Clinical and Laboratory Standards Institute (CLSI). Procedures for Validation of INR and Local Calibration of PT/INR Systems. CLSI document H54-A. Wayne, PA: CLSI (2005)

Clinical and Laboratory Standards Institute (CLSI): Protocol for the Evaluation, Validation and implementation of coagulometers. CLSI document H57-P. Wayne, PA: CLSI (2008)

Clinical and Laboratory Standards Institute (CLSI): Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures. CLSI document EP17-A2. Wayne, PA: CLSI (2012)

Dahmane E, Mercier T, Zanolari B, Cruchon S, Guignard N, Buclin T, Leyvraz S, Zaman K, Csajka C, Decosterd LA. An ultra performance liquid chromatography-tandem MS assay for tamoxifen metabolites profiling in plasma: first evidence of 4'-hydroxylated metabolites in breast cancer patients. *J Chromatogr B Analyt Technol Biomed Life Sci.* 878:3402-14 (2010)

Dati F, Panteghini M, Apple FS, Christenson RH, Mair J, Wu AH. Proposals from the IFCC Committee on Standardization of Markers of Cardiac Damage (C-SMCD): strategies and concepts on standardization of cardiac marker assays. *Scand J Clin Lab Invest Suppl.* (1999)

Dasgupta A, Wahed A. *Clinical Chemistry, Immunology and Laboratory Quality Control. A Comprehensive Review for Board Preparation, Certification and Clinical Practice.* eBook ISBN: 9780124079359 Hardcover ISBN: 9780124078215, Elsevier, 2014. Chapter 2. Immunoassay Platform and Designs 19-34 (Dasgupta A et al, 2014a)

Dasgupta A, Wahed A. *Clinical Chemistry, Immunology and Laboratory Quality Control. A Comprehensive Review for Board Preparation, Certification and Clinical Practice.* eBook ISBN: 9780124079359 Hardcover ISBN: 9780124078215, Elsevier, 2014. Chapter 8. Cardiac markers 127-143 (Dasgupta A et al, 2014b)

Dasgupta A, Wahed A. *Clinical Chemistry, Immunology and Laboratory Quality Control. A Comprehensive Review for Board Preparation, Certification and Clinical Practice.* eBook ISBN: 9780124079359 Hardcover ISBN: 9780124078215, Elsevier, 2014. Chapter 14 Therapeutic Drug Monitoring 249-271 (Dasgupta A et al, 2014c)

Dasgupta A, Wahed A. *Clinical Chemistry, Immunology and Laboratory Quality Control. A Comprehensive Review for Board Preparation, Certification and Clinical Practice.* eBook ISBN: 9780124079359 Hardcover ISBN: 9780124078215, Elsevier, 2014. Chapter 15 Interferences in Therapeutic Drug Monitoring 276-287 (Dasgupta A et al, 2014d)

Decosterd LA, Widmer N, Zaman K, Cardoso E, Buclin T, Csajka C. Therapeutic drug monitoring of targeted anticancer therapies. *Biomark Med.* 9:887-893 (2015)

Derissen EJ, Hillebrand MJ, Rosing H, Schellens JH, Beijnen JH. Development of an LC-MS/MS assay for the quantitative determination of the intracellular 5-fluorouracil nucleotides responsible for the anticancer effect of 5-fluorouracil. *J Pharm Biomed Anal.* 110:58-66 (2015)

Dupouey J, Doudka N, Belo S, Blin O, Guilhaumou R. Simultaneous determination of four antiepileptic drugs in human plasma samples using an ultra-high-performance liquid chromatography tandem mass spectrometry method and its application in therapeutic drug monitoring. *Biomed Chromatogr.* 30:2053-2060 (2016)

Dybkaer R. From total allowable error via metrological traceability to uncertainty of measurement of the unbiased results. *Accred Qual Assur* 4:401-405 (1999)

English E, Lenters-Westra E. HbA1c method performance: The great success story of global standardization. *Crit Rev Clin Lab Sci.* 5:408-419 (2018)

EurA1c Trial Group. EurA1c: The European HbA1c Trial to investigate the performance of HbA1c assays in 2166 laboratories across 17 countries and 24 manufacturers by use of the IFCC model for quality targets. *Clin. Chem.* 64:1183-1192 (2018)

Favaloro EJ, McVicker W, Lay M, Ahuja M, Zhang Y, Hamdam S, Hocker N. Harmonizing the International Normalized Ratio (INR): Standardization of Methods and Use of Novel Strategies to Reduce Interlaboratory Variation and Bias. *Am J Clin Pathol.* 145:191-202 (2016)

- Favarolo EJ. Optimizing the verification of mean normal Prothrombin Time (MNPT) and International Sensitivity Index (ISI) for accurate conversion of Prothrombin Time (PT) to International Normalized Ratio. *Hemostasis and Thrombosis: Methods and Protocols, Methods in Molecular Biology* vol. 1646 DOI 10, 1007/978-1-4939-7198-1\_4. Eds. Favarolo EJ and Lippi G. Springer Science-Business Media LLC 2017
- Favaloro EJ, Gosselin R, Olson J, Jennings I, Lippi G. Recent initiatives in harmonization of hemostasis practice. *Clin Chem Lab Med*. 56:1608-1619 (2018)
- Fraser CG. Reference change values. *Clin Chem Lab Med*. 50: 807-815 (2002)
- Ferraro S, Braga F, Panteghini M. Laboratory medicine in the new healthcare environment. *Clin Chem Lab Med* 54:523-33 (2016)
- Fiaturi N, Greenblatt DJ. Therapeutic Drug Monitoring of Antidepressants. *Handb Exp Pharmacol*. 2018 Sep 8. doi: 10.1007/164\_2018\_161.
- Filler G, Smith N. The need for tacrolimus assay standardization. *Ther. Drug Monit*. 36:693-695 (2014)
- Green SF. The cost of poor blood specimen quality and errors in preanalytical processes. *Clin Biochem* 46: 1175-1179 (2013).
- Font P, Loscertales J, Soto C, Ricard P, Novas CM, Martín-Clavero E, López-Rubio M, Garcia-Alonso L, Callejas M, Bermejo A, Benavente C, Ballesteros M, Cedena T, Calbacho M, Urbina R, Villarrubia J, Gil S, Bellón JM, Díez-Martín JL, Villegas A. Interobserver variance in myelodysplastic syndromes with less than 5 % bonemarrowblasts: unilineage vs. multilineage dysplasia and reproducibility of the threshold of 2 % blasts. *Ann Hematol*. 94:565-73 (2015)
- Fornasaro S, Bonifacio A, Marangon E, Buzzo M, Toffoli G, Rindzevicius T, Schmidt MS, Sergo V. Label-Free Quantification of Anticancer Drug Imatinib in Human Plasma with Surface Enhanced Raman Spectroscopy. *Anal Chem*. 90:12670-12677 (2018)
- Gosselin RC, Adcock DM, Bates SM, Douxfils J, Favaloro EJ, Gouin-Thibault I, Guillermo C, Kawai Y, Lindhoff-Last E, Kitchen S. International Council for Standardization in Haematology (ICSH) Recommendations for Laboratory Measurement of Direct Oral Anticoagulants. *Thromb Haemost*. 118:437-450 (2018)
- Gotta V, Widmer N, Decosterd LA, Chalandon Y, Heim D, Gregor M, Benz R, Leoncini-Francini L, Baerlocher GM, Duchosal MA, Csajka C, Buclin T. Clinical usefulness of therapeutic concentration monitoring for imatinib dosage individualization: results from a randomized controlled trial. *Cancer Chemother Pharmacol*. 74:1307-19 (2014)
- Hanas R, John G. 2010 consensus statement on the worldwide standardization of the haemoglobin A1c measurement. *Diabet. Med*. 27:737-738 (2010)
- Hanas R, John WG; International HbA1c Consensus Committee. 2013 Update on the worldwide standardization of the hemoglobin A1c measurement. *Pediatr Diabetes*. 15:e1-2. (2014)
- Haouala A, Zanolari B, Rochat B, Montemurro M, Zaman K, Duchosal MA, Ris HB, Leyvraz S, Widmer N, Decosterd LA. Therapeutic Drug Monitoring of the new targeted anticancer agents imatinib, nilotinib, dasatinib, sunitinib, sorafenib and lapatinib by LC tandem mass spectrometry. *J Chromatogr B Analyt Technol Biomed Life Sci*. 877:1982-96 (2009)

- Hoelzel W, Weykamp C, Jeppsson JO, Miedema K, Barr JR, Goodall I, Hoshino T, John WG, Kobold U, Little R, Mosca A, Mauri P, Paroni R, Susanto F, Takei I. IFCC reference system for measurement of hemoglobin A1c in human blood and the national standardization schemes in the United States, Japan, and Sweden: a method-comparison study. *Clin. Chem.* 50:166-174 (2004)
- Hoskins JM, Gross AS, Shenfield GM, Rivory LP. High-performance liquid chromatography-electrospray ionization mass spectrometry method for the measurement of moclobemide and two metabolites in plasma. *J Chromatogr B Biomed Sci Appl.* 2001 754:319-26 (2001)
- Hubbard AR, Heath AB. Standardization of factor VIII and von Willebrand factor in plasma calibration of the WHO 5th International Standard (02/150). *Thromb. Hemost.* 2: 1380-1384 (2004)
- Hubbard AR. International biological standards for coagulation factors and inhibitors. *Semin. Thromb. Hemost.* 33: 283-289 (2007)
- Hytest. Single Epitope Sandwich assay for BNP – SES-BNP™. <https://www.hytest.fi/resources/product-info/single-epitope-sandwich-assay-for-bnp-ses-bnp>
- Infusino I, Schumann G, Ceriotti F and Panteghini M. Standardization in clinical enzymology: a challenge for the theory of metrological traceability. *Clin Chem Lab Med* 48:301–307 (2010)
- Infusino I, Fruscicante E, Braga F, Panteghini M. Progress and impact of enzyme measurement standardization. *Clin. Chem. Lab. Med.* 28:155-161 (2016)
- ISO 18153:2003: In vitro diagnostic medical devices – Measurement of quantities in biological samples – Metrological traceability of values for catalytic concentration of enzymes assigned to calibrators and control materials. ISO, Geneva, Switzerland
- ISO/IEC Guide 98:1995. Guide to the expression of uncertainty in measurement (GUM). International Organization for Standardization: Geneva, 1995.
- ISO 15189:2012. Medical laboratories – Requirements for quality and competence. CEN, Brussels, Belgium, 2012.
- Jannetto PJ, Bratanow NC. Utilization of pharmacogenomics and therapeutic drug monitoring for opioid pain management. *Pharmacogenomics.* 10:1157-67 (2009)
- Jiratham-Opas J, Kanjanavanit R, Wongcharoen W, Punyawudho B, Arunmanakul P, Amaritakomol A, Topaiboon P, Gunaparn S, Phrommintikul A. Can available mathematical models predict serum digoxin levels in Thai patients? *J Clin Pharm Ther.* 43:377-384. (2018)
- Kalp Yetersizliği Birliğinin İşbirliğiyle hazırlanmıştır; Heart Failure Association. ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure. *Türk Kardiyol Dern Ars. Suppl* 3:77-137 (2012)
- Kobold U, Jeppsson JO, Duffler T, Finke A, Hoelzel W, Miedema K. Candidate reference methods for haemoglobin A1c based on epitope mapping. *Clin Chem* 43:1944-1951 (1997)
- Lippi G, Guidi GC, Mattiuzzi C, Plebani M. Preanalytical variability: the dark side of the moon in laboratory testing. *Clin Chem Lab Med.* 44:358-365 (2006)
- Little TA. Method validation essentials, limit of blank, limit of detection, and limit of quantitation. *BioPharm International.* 28(4) (2015)

Lopez-Sanchez RDC, Lara-Diaz VJ, Aranda-Gutierrez A, Martinez-Cardona JA, Hernandez HPLC Method for Quantification of Caffeine and Its Three Major Metabolites in Human Plasma Using Fetal Bovine Serum Matrix to Evaluate Prenatal Drug Exposure. *J Anal Methods Chem.* 2018:2085059 (2018)

Lucha S, Taibon J, Pongratz S, Geletneky C, Huber E, Wintterle-Roehm C, Lang R, Grimm SH, Duelffer T, Tarasov K, Zander J, Vogeser M, Kobold U. An LC-MS/MS based candidate reference method for the quantification of total gentamicin in human serum and plasma using NMR characterized calibrator material. *Clin Chim Acta* 464:211-217 (2017)

Luckenbill K.N., Christenson R.H., Jaffe A.S., Mair J, Ordonez-Llanos J., Pagani F., et al., Cross-reactivity of bnp, nt-probnp, and probnp in commercial bnp and nt-probnp assays: preliminary observations from the ifcc committee for standardization of markers of cardiac damage, *Clin. Chem.* 54:619–621 (2008)

Marlar RA, Gausmann JN, Engel J., Validation of hemostasis and coagulation assays: recommendations and guidelines. *Semin. Thromb. Hemost.* 40:186-194 (2014)

Matsuda A, Kawabata H, Tohyama K, Maeda T, Araseki K, Hata T, Suzuki T, Kayano H, Shimbo K, Usuki K, Chiba S, Ishikawa T, Arima N, Nohgawa M, Ohta A, Miyazaki Y, Nakao S, Ozawa K, Arai S, Kurokawa M, Mitani K, Takaori-Kondo A; Japanese National Research Group on Idiopathic Bone Marrow Failure Syndromes. Interobserver concordance of assessments of dysplasia and blastcounts for the diagnosis of patients with cytopenia: From the Japanese central review study. *LeukRes.* 74:137-143 (2018)

McMurray JJ, Adamopoulos S, Anker SD, Auricchio A, Böhm M, Dickstein K, Falk V, Filippatos G, Fonseca C, Gomez-Sanchez MA, Jaarsma T, Kober L, Lip GY, Maggioni AP, Parkhomenko A, Pieske BM, Popescu BA, Ronnevik K, Rutten FH, Schwitter J, Seferovic P, Stepinska J, Trindade PT, Voors AA, Zannad Falez, Zeiher A; Miller WG, Myers GL, Rej R. Why commutability matters. *Clin. Chem.* 52:533-534 (2006)

Menditto A, Patriarca M, Magnusson B. Understanding the meaning of accuracy, trueness and precision. *Accred Qual Assur.* 12:45-47 (2006)

Meroni PL, Biggioggero M, Pierangeli SS, Sheldan J, Zegers I, Borghi MC. Standardization of autoantibody testing: a paradigm for serology in rheumatic diseases. *Nat. Rev. Rheumatol.* 10:35-43 (2014)

Miller WG, Myers GL, Rej R. Why commutability matters. *Clin. Chem.* 52:533-534 (2006)

Miller WG, Schimmel H, Rej R, Greenberg N, Ceriotti F, Burns C, Budd JR, Weykamp C, Delatour V, Nilsson G, MacKenzie F, Panteghini M, Keller T, Camara JE, Zegers I, Vesper HW; IFCC Working Group on Commutability. IFCC Working Group Recommendations for Assessing Commutability Part 1: General Experimental Design. *Clin Chem.* 64:447-454 (2018)

Morishige H, Shuto H, Ieiri I, Otsubo K, Oishi R. Instability of standard calibrators may be involved in overestimating vancomycin concentrations determined by fluorescence polarization immunoassay. *Ther Drug Monit.* 18:80-85 (1996)

Morris RG, Jones TE, Goldsworthy SJ, Wagner TJ, Ho H, Horowitz JD. Suspected DLIS interference in the dimension DGNA digoxin assay method and the clinical application of the revised digoxin target range. *Ther Drug Monit.* 28:454-457 (2006)

Morshed SA, Davies TF, Graves' disease mechanisms: the role of stimulating, blocking, and cleavage region TSH receptor antibodies. *Horm. Metab Res.* 47:880-888 (2015)

Mosca A, Paleari R. Sources and performance criteria of uncertainty of reference measurement procedures. *Clin Biochem.* 57:29-36 (2018)

Nafziger AN, Bertino JS. Utility and application of urine drug testing in chronic pain management with opioids. *Clin J Pain* 25:73-79 (2009)

Neeves KB, Mccarty OJT, Rininger AJ, Sugimoto M, Kino MR, for the Biorheology Subcommittee of the SSC of the ISTH. Flow-dependent thrombin and fibrin generation *in vitro*: opportunities for standardization : communication from SSC of the ISTH. *J. Thromb. Haemost.* 12:418-420 (2013)

Nilsson G, Budd JR, Greenberg N, Delatour V, Rej R, Panteghini M, Ceriotti F, Schimmel H, Weykamp C, Keller T, Camara JE, Burns C, Vesper HW, MacKenzie F, Miller WG; IFCC Working Group on Commutability. IFCC Working Group Recommendations for Assessing Commutability Part 2: Using the Difference in Bias between a Reference Material and Clinical Samples. *Clin Chem.* 64:455-464 (2018)

Nunes LAS, Brenzikofer R, Vaz de Macedo D. Reference change values of blood analytes from physically active subjects. *Eur J Appl Physiol* 110:191-198 (2010)

Olson ID, Brandt JT, Chandler WL, Van Cott EM, Cunningham MT, Hayes TE, Kottke-Marchant KK, Makar RS, Uy AB, Wang EC. Laboratory reporting of the international normalized ratio: progress and problems. *Arch. Path. Lab. Med.* 131:1641-1647 (2007)

Orrico KB, Wu M, Wilson AR. Assessment of the appropriateness of serum digoxin concentration measurement in a medical group setting. *J Manag Care Pharm.* 17:695-700 (2011).

Palmer L, Briggs C, McFadden S, Zini G, Burthem J, Rozenberg G, Proytcheva M, Machin SJ. ICSH recommendations for the standardization of nomenclature and grading of peripheral blood cell morphological features. *Int J Lab Hematol.* 37:287-303 (2015)

Panteghini M. Recent approaches to the standardization of cardiac markers. *Scand. J. Clin. Lab. Invest.* 61:95-102 (2001)

Panteghini M, Ceriotti F, Schumann G, Siekmann L. Establishing a reference system in clinical enzymology. *Clin. Chem. Lab. Med.* 39:795-800 (2001)

Panteghini M. Current concepts in standardization of cardiac marker immunoassays. *Clin Chem Lab Med.* 42:3-8. (2004)

Panteghini M, Ceriotti F, Jones G, Oosterhuis W, Plebani M and Sandberg S, Task Force on Performance Specifications in Laboratory Medicine of the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM). Strategies to define performance specifications in laboratory medicine: 3 years on from the Milan Strategic Conference. *Clin Chem Lab Med.* 55:1849-1856 (2017)

Penttilä I, Penttilä K, Holm P, Laitinen H, Ranta P, Törrönen J, Rauramaa R. Methods, units and quality requirements for the analysis of haemoglobin A1c in diabetes mellitus. *World J Methodol.* 6:133-142 (2016)

Pesce AJ, Rashkin M, Kotagal U. Standards of laboratory practice: theophylline and caffeine monitoring. National Academy of Clinical Biochemistry. *Clin Chem.* 44:1124-1128 (1998)

Plebani M. Errors in clinical laboratories or errors in laboratory medicine? *Clin Chem Lab Med* 44:750-759 (2006)

Plebani M, Astion ML, Barth JH, Chen W, de Oliveira Galoro CA, Escuer MI, et al. Harmonization of quality indicators in laboratory medicine. A preliminary consensus. *Clin Chem Lab Med* 52:951-958 (2014)

- Plebani M. Quality in laboratory medicine: 50 years on. *Clin Biochem* 50:101-104 (2017)
- Plebani M. Analytical quality: an unfinished journey. *Clin Chem Lab Med* 56:357-359.(2018)(Plebani, 2018a)
- Plebani M. Harmonization in laboratory medicine: more than clinical chemistry? *Clin Chem Lab Med* 56:1579-1586 (2018) (Plebani, 2018b)
- Poller L, Ibrahim S, Pattison A, Jespersen J. European Action on Anticoagulation (EAA), formerly European Concerted Action on Anticoagulation (ECAA). INR derivation with the PT/INR Line simplified using a spreadsheet from the world wide web. *J Clin. Pathol.* 64:930-932 (2011)
- Rami L, Roura M, Canalias F. Evaluation of commutability of several materials for harmonization alkaline phosphatase catalytic concentration measurements. *Clin Chim Acta.* 413:1249-1254 (2012)
- Reimers A, Berg JA, Burns ML, Brodtkorb E, Johannessen SI, Johannessen Landmark C' Reference ranges for antiepileptic drugs revisited: a practical approach to establish national guidelines. *Drug Des Devel Ther.* 12:271-280 (2018)
- Rosenberg V, Donald A. Evidence based medicine: an approach to clinical problem-solving. *BMJ* 310:1122-1126 (1995)
- Saenger AK, Rodriguez-Fraga O, Ler R, Ordonez-Llanos J, Jaffe AS, Goetze JP, Apple FS. Specificity of B-Type Natriuretic Peptide Assays: Cross-Reactivity with Different BNP, NT-proBNP, and proBNP Peptides. *Clin Chem.* 63:351-358 (2017)
- Sameri RM, Soberman JE, Finch CK, Self TH. Lower serum digoxin concentrations in heart failure and reassessment of laboratory report forms. *Am J Med Sci.* 324:10-13 (2002).
- Schreiber-Deturmeny E Bruguerolle B. Simultaneous high-performance liquid chromatographic determination of caffeine and theophylline for routine drug monitoring in human plasma. *J Chromatogr B Biomed Appl.* 677:305-312. (1996)
- Schumann G, Bonnora R, Ceriotti F, Clerc-Renaud P, Ferrero CA, Ferard G et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 2. Reference procedures for the measurement of catalytic concentration of creatine kinase. *Clin. Chem. Lab. Med* 40:635-642 (2002) (Schumann G et al, 2002a)
- Schumann G, Bonnora R, Ceriotti F, Clerc-Renaud P, Ferrero CA, Ferard G et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 3. Reference procedures for the measurement of catalytic concentration of lactate dehydrogenase. *Clin. Chem. Lab. Med* 40:643-648 (2002) (Schumann G et al, 2002b)
- Schumann G, Bonnora R, Ceriotti F, Ferard G, Ferrero CA, Franck PF et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 4. Reference procedures for the measurement of catalytic concentration of alanine aminotransferase *Clin. Chem. Lab. Med* 40:718-724 (2002) (Schumann G et al, 2002c)
- Schumann G, Bonnora R, Ceriotti F, Ferard G, Ferrero CA, Franck PF et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5. Reference procedures for the measurement of catalytic concentration of aspartate aminotransferase *Clin. Chem. Lab. Med* 40:725-733 (2002) (Schumann G et al, 2002d)

Schumann G, Bonnora R, Ceriotti F, Ferand G, Ferrero CA, Franck PF et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 6. Reference procedures for the measurement of catalytic concentration of  $\gamma$ -glutamyl transferase. *Clin. Chem. Lab. Med* 40:754-758 (2002) (Schumann G et al, 2002e)

Schumann G, Aok R, Ferrero CA, Ehlers G, Ferand G, Gela FJ et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 8. Reference procedures for the measurement of catalytic concentration of  $\gamma$ -glutamyl transferase. *Clin. Chem. Lab. Med* 44:1146-1155 (2006)

Schumann G, Ceriotti F, Panteghini M. Standardization in clinical enzymology: a challenge for the theory of metrological traceability. *Clin. Chem. Lab. Med.* 48:301-307 (2010)

Schumann G, Klauke R, Canalias F, Bossert –Reuther S, Franck PF, Gella FJ et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 9. Reference procedures for the measurement of catalytic concentration of alkaline phosphatase *Clin. Chem. Lab. Med* 49:1439-1446 (2011)

Sciacovelli L, Panteghini M, Lippi G, Sumarac Z, et al. Defining a roadmap for harmonizing quality indicators in Laboratory Medicine: a consensus statement on behalf of the IFCC Working Group „Laboratory Error and Patient Safety” and EFLM Task and Finish Group „Performance specifications for the extra-analytical phases”. *Clin Chem Lab Med* 55:1478-1488 (2017)

Seeger C, Shipkova M, Christians U, Billaud EM, Wang P, Holt DW, Brunet M, Kunicki PK, Pawiński T, Langman LJ, Marquet P, Oellerich M, Wieland E, Wallemacq P. Assuring the Proper Analytical Performance of Measurement Procedures for Immunosuppressive Drug Concentrations in Clinical Practice: Recommendations of the International Association of Therapeutic Drug Monitoring and Clinical Toxicology Immunosuppressive Drug Scientific Committee. *Ther Drug Monit.* 38:170-189 (2016)

Semenov AG, Tamm NN, Apple FS, Schulz KM, Love SA, Ler R, Feygina EE, Katrukha AG. Searching for a BNP standard glycosylated proBNP as a common calibrator enables improved comparability of commercial BNP immunoassays. *Clin. Biochem.* 50:181-185 (2017)

Senent L, Arenillas L, Luño E, Ruiz JC, Sanz G, Florensa L. Reproducibility of the World Health Organization 2008 criteria for myelodysplastic syndromes. *Haematologica.* 98:568–575 (2013)

Simons FE, Rigatto H, Simons KJ. Pharmacokinetics of theophylline in neonates. *Semin Perinatol.* 5:337-345 (1981)

Simons SA, Molinelli AR, Sobhani K, Ramey PM, Eloofnagle A. Two cases with unusual vancomycin measurements. *Clin. Chem. Acta* 55:578-582 (2009)

Snell NJC. Drug interactions with anti-asthma medication. *Respir Med* 88:83-88 (1994)

Solberg HE. The IFCC recommendations on estimation of reference intervals. The RefVal Program. *Clin Chem Lab Med* . 42:710-714 (2004)

Soni K, Patel N, Singh K, Jha A, Patel H, Gupta R, Srinivas NR. A Sensitive Triple Quadrupole Liquid Chromatography Mass Spectrometric Method for the Estimation of Valproic Acid in K2EDTA Human Plasma using Furosemide as the Internal Standard. *Drug Res (Stuttg).* 66:666-672. (2016)

Stenman UH. Immunoassay Standardization: Is It Possible, Who Is Responsible, Who Is Capable? *Clin. Chem* 47:815-820 (2001)

- Sturgeon C. Standardization of tumor markers – priorities identified through external quality assesment. *Scand. J. Clin. Lab. Invest.* 245:S94-99 (2016)
- Stöckl D, Franzini C, Kratochvila J, Middle J, Ricos C, Siekmann L, Thienpont LM. Analytical specifications of reference method compliation and critical discussion (from the members of the European EQA-Organizers Working Group B) *Eur J Clin Chem Clin Biochem* 34:319-337 (1996)
- Taibon J, Schmid R, Lucha S, Pongratz S, Tarasov K, Seger C, Timm C, Thiele R, Herlan JM, Kobold U. An LC-MS/MS based candidate reference method for the quantification of carbamazepine in human serum. *Clin Chim Acta.* 472:35-40 (2017)
- Tamm NN, et al. Novel immunoassay for quantification of brain natriuretic peptide and its precursor in human blood. *Clin. Chem.* 54:1511-1518 (2008)
- Tang W.H., Francis G.S., Morrow D.A., Newby L.K., Cannon C.P., Jesse R.L., et al., Nationalacademy of clinical biochemistry laboratory medicine practice guidelines: clinical utilization of cardiac biomarker testing in heart failure, *Circulation* 116 e99–109 (2007)
- Thappali SR, Varanasi KV, Veeraraghavan S, Vakkalanka SK, Khagga M. Simultaneous determination of methotrexate, dasatinib and its active metabolite N- deshydroxyethyl dasatinib in rat plasma by LC-MS/MS: method validation and application to pharmacokinetic study. *Arzneimittelforschung.* 62:624-630 (2012)
- Tate J, Panteghini M. Standardization of cardiac troponin I measurements – the way forward. *biochimica clinica,* 35:461-464 (2011)
- Thomas JB, Yen JH, Schantz MM, Porter BJ, Sharpless KE. Determination of caffeine, theobromine, and theophylline in standard reference material 2384, baking chocolate, using reversed-phase liquidchromatography. *J Agric Food Chem.* 52:3259-3263 (2004)
- Thomas JB, Sharpless KE, Mitvalsky S, Roman M, Yen J, Satterfield MB. Determination of caffeine and caffeine-related metabolites in ephedra-containing standard reference materials using liquidchromatography with absorbance detection and tandem mass spectrometry. *J AOAC Int.* 90:934-940 (2007)
- Thomson JM, Taberner DA, Poller L. Automation and prothrombin time: a United Kingdom field study of two widely used coagulometers. *J. Clin. Pathol.* 43:679-684 (1990)
- Tong Q, Chen B, Zhang R, Zuo C. Standardization of clinical enzyme analysis using frozen human serum pools with values assigned by the International Federation of Clinical Chemistry and Laboratory Medicine reference measurement procedures. *Scand J Clin Lab Invest.* 78:74-80 (2018)
- Tonutti E, Visentini D, Bizzaro N. Interpretative comments on autoantibody tests. *Autoimmun Rev.* 6:341-346 (2007)
- Torma AF, Groves K, Biesenbruch S, Mussell C, Reid A, Ellison S, Cramer R, Quaglia M. A candidate liquid chromatography mass spectrometry reference method for the quantification of the cardiac marker 1-32 B-type natriuretic peptide. *Clin Chem Lab Med.* 55:1397-1406 (2017)
- Tripodi A., Chantarangkul V, Pengo V. Standardization of lupus anticoagulant. The Lupus Anticoagulant Sensitivity Index (LASI). *Lupus* 21:715-717 (2012)
- Tricyclic antidepressants – Blood level measurements and clinical outcome: An APA Task Force Report. Task Force on the use of laboratory methods in Psychiatry. *Psychiatry* 142:155-162 (2006)

- Tuzzoli R, Vilalta D, Bizzaro N. Challenges in the standardization of autoantibody testing: a comprehensive review. *Clin Rev Allerg. Immunol.* 53: 68-77 (2017)
- Tozzoli R, Villalta D, Bizzaro N. Challenges in the Standardization of Autoantibody Testing: a Comprehensive Review. *Clin Rev Allergy Immunol.* 53:68-77 (2017)
- Twernbold R, Reichlin T, Reiter M, Muller C. High-sensitive cardiac troponin: friend or foe? *Swiss. Med. Wkly* 141:w3202 (2011)
- Yoong W, Yeo CP. Clinical laboratory diagnostics of immunosuppressants: One laboratory's journey. *Proceedings of Singapore Healthcare* 24:219-224 (2015)
- Van Eeckhaut A, Lanckmans K, Sarre S, Smolders I, Michotte Y, Vesper HW, Myers GL, Miller WG. Validation of bioanalytical LC-MS/MS assays: evaluation of matrix effects. *J Chromatogr B Analyt Technol Biomed Life Sci.* 877:2198-207 (2009)
- Vásárhelyi B. Hemoglobin-A-1c-szint-mérés: analitikai vonatkozások és ezek jelentősége a klinikai döntéshozatalban. *Orv. Hetil.* 157:753-757 (2016)
- Vesper HW, Myers GL, Miller WG. Current practices and challenges in the standardization and harmonization of clinical laboratory tests. *Am J Clin Nutr* 104:907S–12S (2016)
- Vogesser M, Seger C. Pitfalls associated with the use of liquid chromatography-tandem mass spectrometry in the clinical laboratory. *Clin. Chem.* 56:1234-1244 (2010)
- Vu RL, Helmeste D, Albers L, Reist C. Rapid determination of venlafaxine and O-desmethylvenlafaxine in human plasma by high-performance liquid chromatography with fluorimetric detection. *J Chromatogr B Biomed Sci Appl.* 703:195-201 (1997)
- Walz B, Fierz W. The concept of reference change values (RCV). Will it supersede reference intervals? *Ther Umsch.* 72:130-135 (2015)
- Westgard, J.O. Barry, PL, Hunt MR, Groth, T. A Multi-Rule Shewhart Chart For Quality Control In Clinical Chemistry. *Clin. Chem.* 27:493-501 (1981).
- Westgard JO, Westgard SA. Measuring analytical quality. Total analytical error versus measurement uncertainty. *Clin Lab Med* 37:1-13 (2017)
- Weykamp CW, Penders TJ, Muskiet FA, van der Slik W. Effect of calibration on dispersion of glycohemoglobin values determined by 111 laboratories using 21 methods. *Clin. Chem.* 40:138-144 (1994)
- Widmer N, Bardin C, Chatelut E, Paci A, Beinen J, Leveque D, Veal G, Astier A. Review of therapeutic drug monitoring of anticancer drugs part two – Targeted therapies. *Eur J Cancer* 50:2020-2036 (2014)
- Wilson D, Johnston F, Holt D, Moreton M, Engelmayer J, Gaulier J-M, Luthe H, Marquet P, Moscato D, Oellerich M, Mosso R, Streit F, Brunet M, Filler C, Schmid R, Wallemacq P, Barnes G. Multi-center evaluation of analytical performance of the microparticle enzyme immunoassay for sirolimus. *Clin. Biochem.* 39:378-389 (2006)
- Winter WE, Flax SD, Harris NS. Coagulation Testing in the Core Laboratory. *Lab Med.* 48:295-313 (2017)

Whicher JT, Ritchie RF, Johnson AM, Bauchner S, Benvenu J, Blirup-Jensen S, Carlstrom A, Dati F, Ward AM, Svendsen PJ. New international reference preparation for proteins in human serum (RPPHS). *Clin. Chem.* 40:934-938 (1994)

White S, Wong SHY. Standards of laboratory practice: analgesic drug monitoring. *Clinical Chemistry* 44:1110-1123 (1998)

Woods AI, Sanchez-Luceros A, Bermejo E et al. Identification of pW2461, as a novel mutation in the gene responsible for platelet-type von Willebrand disease. *Semin. Thromb. Hemost.* 40:151-160 (2014)

Wu JT, Pieper JK, Wu LH, Peters JL, Isolation and characterization of myoglobin and two major isoforms from sheep heart. *Clin. Chem.* 35:778-782 (1989)

Zhang Y, Zhang R. Recent advances in analytical methods for the therapeutic drug monitoring of immunosuppressive drugs. *Drug Test. Anal.* 10:81-94 (2018)

Zhao L, Yang P, Li P, Wang X, Qin W, Zhang X. Efficiency of individual dosage of digoxin with calculated concentration. *Clin Interv Aging.* 9:1205-1210. (2014)

Zheng B, Li E, Zhu H, Lu J, Shi X, Zhang J, Li M. Automated antinuclear immunofluorescence antibody analysis is a reliable approach in routine clinical laboratories. *Clin Chem Lab Med.* 55:1922-1930. (2017)